



FLOOR BALANCE

LFRB-A20

Floor Balance LFRB-A20 is an electronic balance comprising indicator and square weighing pan, offers 50 kg maximum weighing capacity. Features 304 stainless steel iron pan with ABS house and aluminum sensor, detachable indicator provides kg, g, oz, lb as weighing units. Equipped with LCD screen display and RS232 interface, this equipment offers counting, percentage, full range tare, animal weighing function with over-load protection, zero tracking and low-battery indicator. With quick response, accurate weighing and high quality, this equipment brings great convenient to user.

Features :

- ❑ Electronic balance with detachable indicator and square weighing pan
- ❑ 304 stainless steel iron pan with ABS house and aluminum sensor
- ❑ Weighing Units: kg, g, oz, lb
- ❑ LCD screen display with white backlight and RS232 interface
- ❑ Offers counting, percentage, full range tare, animal weighing function
- ❑ Over-load protection, zero tracking and low-battery indicator
- ❑ Quick response, accurate weighing and high quality
- ❑ Memory of last setting for convenience

Application :

Floor electronic balances are used in the food, pharmaceutical, textile, geological and jewellery industries, as well as in schools, hospitals, laboratories, medical, chemical and other large factories.

Specifications:

Model	LFRB-A20
Capacity	50 kg
Accuracy	1 g
Speed	1 s

Value stable time	2 s
Display	LCD screen display
Pan type	Square
Weighing units	kg, g, oz, lb
Functions	Counting, percentage, full range tare, animal weighing function, zero tracking
Overload protection	Yes
Indicator	Low battery
Sensor	Aluminum sensor
Pan to floor height	120 mm
Material	304 Stainless steel Iron, ABS house
Power	AC
Pan dimensions (LxB)	400x500 mm
Indicator dimensions (LxH)	250x150 mm
Packaging dimensions (WxDxH)	490x450x170 mm
Net weight	6 kg
Gross weight	7 kg

Standard Accessories:

Accessories no.	Accessories Name
1	AC adapter
2	Rechargeable battery