

Manual Dual Door
Combustible

Cabinets
LCBC-A Series



114 L Combustible Cabinet LCBC-A10

114 L Combustible Cabinet LCBC-A10 is a dual walled manual fire resistance cabinet with 38 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features one adjustable shelf to accommodate storage for different containers with 2 inches vents integrated with flame arresters. Equipped with spill catcher shelves for preventing incidental drips. Inner and outer lead free powder coating provides durability and chemical resistance preventing corrosion and humidity.

Features

- ▶ Capacity - 114 L
- ▶ Dual walled manual
- ▶ 38 mm insulating air space
- ▶ 180° fully opening door
- ▶ Three point latch and manual lock
- ▶ One adjustable shelf
- ▶ Adjustable spill retaining shelves

170 L Combustible Cabinet LCBC-A11

170 L Combustible Cabinet LCBC-A11 is a dual walled manual fire resistance cabinet with 38 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features two adjustable shelves to accommodate storage for different containers with 2 inches vents integrated with flame arresters. Equipped with spill catcher shelves for preventing incidental drips. Inner and outer lead free powder coating provides durability and chemical resistance preventing corrosion and humidity.

Features

- ▶ Capacity - 170 L
- ▶ Dual walled manual
- ▶ 38 mm insulating air space
- ▶ 180° fully opening door
- ▶ Three point latch and manual lock
- ▶ Two adjustable shelves
- ▶ Adjustable spill retaining shelves

Applications

Used for storage of flammable and combustible liquids away from incompatible materials (e.g., oxidizers)

Specification

Model No	LCBC-A10	LCBC-A11
Capacity	114 L	170 L
Door type	Dual door manual	Dual door manual
Adjustable shelf	1	2
Shelves loading capacity	50 kgs	100 kgs
Dimension (H x W x D)	1120 x 1090 x 460 mm	1650 x 1090 x 460 mm
Weight	107 kgs	148 kgs

227 L Combustible Cabinet LCBC-A12

227 L Combustible Cabinet LCBC-A12 is a dual walled manual fire resistance cabinet with 38 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features two adjustable shelves to accommodate storage for different containers with 2 inches vents integrated with flame arresters. Equipped with spill catcher shelves for preventing incidental drips. Inner and outer lead free powder coating provides durability and chemical resistance preventing corrosion and humidity.

Features

- Capacity - 227 L
- Dual walled manual
- 38 mm insulating air space
- 180° fully opening door
- Three point latch and manual lock
- Two adjustable shelves
- Adjustable spill retaining shelves

340 L Combustible Cabinet LCBC-A12

340 L Combustible Cabinet LCBC-A12 is a dual walled manual fire resistance cabinet with 38 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features two adjustable shelves to accommodate storage for different containers with 2 inches vents integrated with flame arresters. Equipped with spill catcher shelves for preventing incidental drips. Inner and outer lead free powder coating provides durability and chemical resistance preventing corrosion and humidity.

Features

- ▶ Capacity - 340 L
- ▶ Dual walled manual
- ▶ 38 mm insulating air space
- ▶ 180° fully opening door
- ▶ Three point latch and manual lock
- ▶ Two adjustable shelves
- ▶ Adjustable spill retaining shelves

Applications

Used for storage of flammable and combustible liquids away from incompatible materials (e.g., oxidizers)

Specification

Model No	LCBC-A12	LCBC-A13
Capacity	227 L	340 L
Door type	Dual door manual	Dual door manual
Adjustable shelf	2	2
Shelves loading capacity	100 kgs	100 kgs
Dimension (H x W x D)	1650 x 860 x 860 mm	1650 x 1090 x 860 mm
Weight	173 kgs	186 kgs