

Manual Dual Door
Strong Acid and
Alkali Chemical
Storage Cabinet
LSAC-A1 SERIES



Manual Dual Door Strong Acid and Alkali Chemical Storage Cabinet

105 L Strong Acid and Alkali Cabinet LSAC-A10



105 L Strong Acid and Alkali Cabinet LSAC-A10 is a dual walled acid and alkali cabinet with 50 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features one adjustable shelf to accommodate storage for different containers. Inner and outer epoxy resin powder coating provides durability and chemical resistance preventing corrosion and humidity.

Features

- Capacity - 105 L
- Dual walled manual
- 50 mm insulating air space
- Three point latch and manual lock
- One adjustable shelf

Applications

Used to store all kinds of corrosive material or weak acid and alkali, for example: acetic acid, Carbonic acid, phosphoric acid, silicic acid, etc.

Specifications

Model No	LSAC-A10
Capacity	105 L
Door type	Dual door manual
Adjustable shelf	1
Dimension (W x D x H)	910 x 600 x 900 mm
Weight	60 kgs

Manual Dual Door Strong Acid and Alkali Chemical Storage Cabinet

114 L Strong Acid and Alkali Cabinet LSAC-A11



114 L Strong Acid and Alkali Cabinet LSAC-A11 is a dual walled acid and alkali cabinet with 50 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features one adjustable shelf to accommodate storage for different containers. Inner and outer epoxy resin powder coating provides durability and chemical resistance preventing corrosion and humidity.

Features

- Capacity - 105 L
- Dual walled manual
- 50 mm insulating air space
- Three point latch and manual lock
- One adjustable shelf

Applications

Used to store all kinds of corrosive material or weak acid and alkali, for example: acetic acid, Carbonic acid, phosphoric acid, silicic acid, etc

Specifications

Model No	LSAC-A11
Capacity	114 L
Door type	Dual door manual
Adjustable shelf	1
Dimension (W x D x H)	1090 x 460 x 1120 mm
Weight	70 kgs

Manual Dual Door Strong Acid and Alkali Chemical Storage Cabinet

170 L Strong Acid and Alkali Cabinet LSAC-A12



170 L Strong Acid and Alkali Cabinet LSAC-A12 is a dual walled acid and alkali cabinet with 50 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features one adjustable shelf to accommodate storage for different containers. Inner and outer epoxy resin powder coating provides durability and chemical resistance preventing corrosion and humidity.

Features

- Capacity - 170 L
- Dual walled manual
- 50 mm insulating air space
- Three point latch and manual lock
- One adjustable shelf

Applications

Used to store all kinds of corrosive material or weak acid and alkali, for example: acetic acid, Carbonic acid, phosphoric acid, silicic acid, etc.

Specifications

Model No	LSAC-A12
Capacity	170 L
Door type	Dual door manual
Adjustable shelf	1
Dimension (W x D x H)	1090 x 460 x 1650 mm
Weight	102 kgs