



MEDICAL REFRIGERATOR

LMR-B15

MEDICAL REFRIGERATOR LMR-B15

Medical Refrigerator LMR-B15 is upright unit microprocessor temperature controlled unit, refrigerator with user-friendly control panel and a LED display for clear observation of temperature readings. Designed with 130 L capacity, 2 to 8°C temperature range, R600a CFC free refrigerant, and a forced air refrigeration cooling circulation system. Features auto defrost mode, adjustable rust proof shelves, single glass door with electric heating. Offers audible and visual alarm for high and low temperature, Power failure, Sensor failure, Door ajar

FEATURES

- LED digital display of inner temperature and ambient temperature
- Managed and monitored by an integrated microprocessor controller
- Single glass door design, safety lock with key for safe operation
- Easy to clean, corrosion resistant, stainless steel lining
- Built in with multiple visual and audible alarm functions
- Forced air refrigeration cooling system
- Four pieces of Adjustable shelves
- Four moveable casters with lock for easy mobility and handling
- Auto defrost mode
- Audible and Visual alarm for high and low temperature, Power failure, Sensor failure, Door ajar

APPLICATIONS

Medical Refrigerator used in storage of drugs, vaccines, reagents and biological products. It is suitable to use in hospitals, clinics, research institutions, universities and colleges, military enterprises, laboratories, centres for disease prevention & control, pharmaceutical factories.

SPECIFICATION

Model no.	LMR-B15
Capacity	130 L
Temperature Rang	2 to 8°C
Temperature Accuracy	0.1°C
Refrigerant	R600a, CFC free
Refrigeration Type	Forced Air Refrigeration Cooling System
Condenser And Evaporator	Bundy Tube Condenser
Defrost	Auto Defrost Mode
Display	LED
Control System	Microprocessor Temperature Controlled
Shelves	04 Pieces Adjustable
Door Type	Single Glass Door With Electric Heating
Power Consumption	136 W
Power Supply	AC1108220V ±10%, 50/60 Hz
External Dimension (W×D×H)	500 × 550 × 1300 mm
Packaging Dimension (W×D×H)	590 × 595 × 1380 mm
Gross weight	51 kg

STANDARD ACCESSORIES

- LED lamp
- 4 Shelves
- Door Keys
- USB port

