

 LABTRON



# Solar Freezers

## ■ Solar Freezer LSF-A Series

---

These chest-type Solar Freezer are eco-friendly, safe, portable, high-performing, fast cooling and high temperature resistant, with advance design, low maintenance for quality performance, safety and accurate results.



*LSF-A10*

### ■ FEATURES:

---

- Environment friendly CFC free refrigerant R134a is used
- Design with single top door with door lock to prevent unauthorized access
- Automatic low voltage disconnection for protection
- Low maintenance and easy to clean
- Fast cooling with optional temperature control
- Freezer is equipped with drawers for safe storage
- Designed to ensure high temperature resistance and minimum energy consumption
- Equipped with hermetic domus compressor with integrated electronic control panel

## APPLICATION:

Widely used in body fluid banks, food industries, life science, biotechnology, genetics, research and development.

## TECHNICAL PARAMETERS:

Model	LSF-A10	LSF-A11	LSF-A12	LSF-A13	LSF-A14
Capacity	100Litres	150Litres	200Litres	268Litres	358Litres
Input Voltage	DC12V/24V				
Power Consumption	65W	75W	87W	100W	120W
Product Size (L*W*H mm)	580*610*870	580*610*870	910*610*870	1100*610*870	1344*610*870
Carton Size (L*W*H mm)	580*610*870	580*610*870	910*610*870	1100*610*870	1344*610*870
N.W./G.W(Kg)	30/33	35/38	42/45	48/51	48/51
Additional Parts	Solar panel: 200W; Battery: 150Ah; Charge controller: 12V/24V-15A				



Labtron Equipment Ltd

Basing View, Basingstoke, Hampshire RG21 4RG, UK

+44 1256 610 110 | [info@Labtron.com](mailto:info@Labtron.com) | [www.Labtron.com](http://www.Labtron.com)

